

## Release Notes LogicNets AUC Solution Version October 2023

---

### Contents

1.0	Overview .....	1
2.0	LogicNets AUC Version Number .....	1
3.0	AUC Solution Enhancements .....	1
3.1	Updated LogicNets Release .....	1
3.2	AUC Validation Reconfiguration.....	1
3.3	Additional ICD-10 Mappings .....	1

### 1.0 Overview

LogicNets has released a new version of the LogicNets AUC Solution. In addition to migrating the AUC Solution to the latest version of the LogicNets application, LogicNets updated elements of the AUC Solution to streamline the user experience.

### 2.0 LogicNets AUC Version Number

The new release is version number 3.0.35988 (8.0.0.35988).

### 3.0 AUC Solution Enhancements

LogicNets made the following updates to the AUC Solution.

#### 3.1 Updated LogicNets Release

LogicNets updated the AUC Solution to run on our latest version release, Version 8.0. This version includes performance and security improvements resulting from upgrades to LogicNets' underlying software infrastructure.

#### 3.2 AUC Validation Reconfiguration

With this release, LogicNets reconfigured the AUC validation service to run in "Relaxed Mode" as the default. This allows the system to return a validation response of "appropriate" for a single advanced imaging service, even though multiple imaging services may be recommended to meet the clinical guidance. For example, an MRI and CT may be required to meet full recommended guidance, but Relaxed Mode allows ordering professionals to get validation for the CT order alone.

#### 3.3 Additional ICD-10 Mappings

This release of the AUC Solution includes additional ICD-10 code input mappings for guideline selection filtering. The guideline filtering feature to which this applies is generally activated by administrators for organizations accessing both the NCCN and Intermountain Health (IHS) guidelines.